

m.

0

.

Camtin AR POLLUTION CONTROL

Buld

GOLD SERIES CAMTAIN

DUST COLLECTOR FOR PHARMACEUTICAL APPLICATIONS

camfilapc.com

CAMFIL APC – EXPERTS FOR PHARMA APPLICATIONS

Camfil Air Pollution Control is part of Camfil, the world's largest air filtration manufacturer, with presence in more than 150 countries and production on all continents. We are experienced in pharmaceuitcal support, with a team of professionals focused solely on pharmaceutical applications, making us your ideal partner. We work together with pharmaceutical and engineering companies and contractors to find effective, efficient, reliable solutions to control dust in pharmaceutical operations to protect workers, the environment and equipment while maximizing production.

GOLD SERIES CAMTAIN

Camfil APC has named its contained dust collection system the Gold Series Camtain. Camfil is the parent company of Camfil Air Pollution Control and is the world

Safe-change containment systems are available for both the filter cartridges and discharge system. The cartridge

change utilizes a bag-in/bag-out (BIBO) method while

Gold Series can also support traditional dust collection

for nuisance dusts and fumes that do not require full

isolation and containment. Camfil APC has contained

and noncontained dust collectors for pharmaceutical

applications in the Americas, Asia and Europe.

thedischarge uses continuous liner technology. The

leader in contained HEPA filtration systems.



Gold Series Camtain



Potent Compound Containment Collectors

BENEFITS OF THE GOLD SERIES COLLECTOR

- ✓ High entry, cross flow inlet eliminates upward velocities that can hold fine particulates in the filter cartridges, reducing the re-entrainment of the particulate matter.
- ✓ Vertically arranged filters shed virtually all the particles versus horizontal filters which allow the particles to build on top of the filter.
- ✓ High efficiency filters stop 99.99% of the dust at 0.5 microns!
- ✓ Specially treated filter media repels fine particulates for a lower pressure drop and longer filter life.
- ✓ HemiPleat Gold Cone design provides 25% more media for long service life
- ✓ Filter cartridges are sealed via an internal cam lock action, allowing convenient change-out through the BIBO safechange system.

GOLD SERIES FOR PHARMACEUTICAL MANUFACTURING APPLICATIONS

The Gold Series (GS) units can be used in a variety of pharmaceutical dust collection applications including tablet presses, coating, fluid bed drying, spray drying, blending, granulation and general room ventilation. The GS is perfect for high efficiency filtration for pharmaceutical manufacturing.

Serie

Camfi

.

SURROGATE TESTED DUST COLLECTION SYSTEM FOR PERFORMANCE VERIFICATION

The Gold Series Camtain contained dust collection system has been independently surrogate tested for validated performance verification. A major pharmaceutical company's surrogate testing protocol was followed with an independently contracted, AIHA accredited laboratory (Bureau Veritas) performing the testing. Using milled lactose as the surrogate, we collected over 48 area, personal and swab samples for both the BIBO cartridge filter change and the continuous liner discharge. The GS Camtain can handle risk based category 3, 4 and 5 compounds with OELs less than 1.0 mcg/m3 for an 8-hour time weighted average. Full test report data is available upon request.





HEMIPLEAT® FILTERS WILL IMPROVE THE PERFORMANCE OF ANY CARTRIDGE DUST COLLECTOR.

Power up your dust collector with HemiPleat. With lower pressure drop HemiPleat filters, you can pull more air with less energy, thus capturing pollutants better. Filtration efficiencies exceed 99.99% at 0.5 micron particle by weight.



GOLD CONE TECHNOLOGY

The patented Gold Cone filter has allowed many facilities to reduce the number of filters they have to use and change. The innovative cone of filter media expands the usable area of the filter, reducing the required number of filters by at least a third. The design also promotes long filter life with low pressure drop.



THE GOLD SERIES CAMTAIN DUST COLLECTOR COMBINES ENHANCED PERFORMANCE, SAFE-CHANGE CAPABILITY AND EASE OF MAINTENANCE WHILE PROTECTING THE WORKPLACE AND ENVIRONMENT FROM HARMFUL DUSTS.

Two key concerns are the focus when handling pharmaceutical dusts - the potent, toxic or allergenic properties of the compound as it relates to personnel exposure and the explosion properties of the compound.

The first issue involves understanding the toxicological properties of the material, reviewing the Occupational Exposure Limit (OEL) and performing a risk based exposure evaluation to determine the methods for proper control. In most cases, some level of isolation and containment is required due to the fact that the pharmaceutical dust is extremely potent while being captured in a non-production area and cannot be released into the surrounding environment. In most cases, Camfil APC recommends a HEPA secondary safety system. With HEPA systems after the dust collector, recirculation of the filtered air back into the HVAC system is an option. This can significantly reduce energy costs while providing the necessary level of filtration for discharge air required by OSHA/EPA.

The second concern involves deflagration and explosion potential. Control measures such as explosion venting, chemical suppression and isolation systems may be required depending on the physical characteristics of the dust relating to Kst, Minimum Ignition Energy (MIE) and the location of the collector. When explosion vents are required, they must be vented to the outside by either placing the collector outdoors or ducting the vent exhaust to a safe location at a specified distance through the building structure. Camfil APC recommends an independent PE specify what explosion protection is required for a given material as it relates to standards in NFPA, ATE X and the insurance carriers.

amfi

camfil

camfilapc.com

Camfil APC GmbH | Ludwigstaler Str. 149 | 78532 Tuttlingen / Germany Phone: +49 7461 7011-0 | Fax: +49 7461 7011-133 | E-Mail: apcde-info@camfil.com