



## Advantages

- Designed specifically for pharmaceutical and containment applications
- High collector efficiency using HemiPleat cartridges
- Customised for Original Equipment Manufacturers (OEM)
- Simple cartridge replacement using quick release cam bars
- Bag-in/bag-out safe change options available.
- Modular design for optimum flexibility
- Easy to install and maintain
- Up to 25% smaller

**Application:** The Farr Gold Series® Camtain® is used in a wide range of pharmaceutical applications including tablet presses, coating, fluid bed and spray drying, blending, granulation and general ventilation. Contact Camfil for more information

**Type:** Dust collector

**Options:** A wide variety of options are available including: BIBO (bag in-bag out) for Pharmaceutical Applications, Explosion Venting, Special Inlet Designs, Custom Colours, Stainless Steel Construction, Alternative Hopper Designs etc. please contact us with your specific requirements

**Cartridges:** Vertically mounted to shed dust readily for efficient cleaning and longer service life. High filtration efficiency meeting the 5 mg/m<sup>3</sup> or less emissions required to re-circulate the air back into the work place on non hazardous dusts

Safe-change containment systems are available for both the filter cartridges and discharge system underneath the collector.

The cartridge change utilizes the safe change filter replacement method while the discharge uses continuous liner technology.

The Farr Gold Series Camtain is perfect for high efficiency filtration in pharmaceutical manufacturing processes where recovery of the product is not required.

The only dust collector that is potent compound surrogate tested for validated performance verification. Test report available upon request.

