IMPREGNATED CARBON FOR ETHYLENE OXIDE





ADVANTAGES

- ISO 10121-1:2014 tested media
- Available in a range of sizes
- Activated carbon impregnated to target ethylene oxide
- Able to be installed in any manufacturer's loose-fill filters

Application	Targets ethylene oxide. Used mostly for sterilization in the biopharmaceutical industry.
Media	Impregnated Activated Carbon
Max Temperature (°C)	40
Min Temperature (°C)	10
Relative Humidity max	40%-90%
Installation Options	All loose-fill molecular air filters.

